

CON30

Conductivity Tester

(Conductivity / TDS / Salinity)

2014
New
Model



Precision easy operation: all functions operated in one hand.

Large LCD with backlight, Multi-line display for easy reading,
Auto-Power Off saves battery after 5 mins non-use.

1*1.5V AAA long battery life

Easy maintenance, simple to change battery and electrode

Floats on water, IP67 rated

Field throw measurement (Auto Lock Function)



CLEAN CON30 Conductivity Tester

- Broad measuring range: 0.0 $\mu\text{S}/\text{cm}$ - 20.00 mS/cm , from pure water to waste water applications.
- 4 parameters in one meter, simple calibration.

Range	0.0 $\mu\text{S}/\text{cm}$ (ppm) - 20.00 mS/cm (ppt)
Resolution	0.1 $\mu\text{S}/\text{cm}$ (ppm) - 0.01 mS/cm (ppt)
Accuracy	$\pm 2\%$ F.S.
Temp. Range	0 - 100.0 $^{\circ}\text{C}$ / 32 - 212 $^{\circ}\text{F}$
Working Temperature	0 - 60.0 $^{\circ}\text{C}$ / 32 - 140 $^{\circ}\text{F}$
Auto Temp. Compensation	0 - 60 $^{\circ}\text{C}$
Calibration	1 or 2 point auto calibrate (0% zero oxygen or 100% in air)
TDS Range	0.0 mg/L (ppm) - 20.00 g/L (ppt)
TDS Coefficient	0.40 - 1.00 Adjustable (Coefficient: 0.5)
Salinity Range	0.0 mg/L (ppm) - 13.00 g/L (ppt)
Salinity Coefficient	0.65
Atmospheric	
Pressure Compensation	600 - 1100 mbar
Auto backlight off	1 minutes
Auto power off	5 minutes without key be pressed
Power	1x1.5V AAA battery
Dimensions	(H×W×D) 185×40×48 mm
Weight	95g
Protection	IP67

CON30 is the best value salt meter in your pocket, with 3 parameters (Conductivity / TDS / Salinity tester) in one and easy calibration. This dip-style tester is ideal for spot checking with greater flexibility.

Package Sets



www.cleaninst.com
Water Analysis Solutions