

TL220 is powerful versatile and indispensable centrifuge for routine laboratories. It is not only ideal for small laboratory; it can be also very useful in medium and large laboratories. It is high performance motor, **Maintenance free, Brushless and Microprocessor** for preselecting of speed and time.

- Easy entry of key pad
- Easy to read, digital display
- Display of correct speed & time volume
- Impulse key for short centrifuge
- Running time in minute, max 99 min
- Open lid key
- Lid locking
- Motor over heating protection



MODEL	TH 220		TH200	
Max Speed	16000r/min		16000 r/min	
Max RCF	23669Xg		17800xg	
Timer Range	0~99min		0~99min	
Power Supply	AC220V50Hz5A		AC220V50Hz5A	
Dimensions(LxWxH)	330x420x275mm		230x310x240mm	
Weight	20kg		15kg	
	Max Speed (r/min)	MaxRCF ( xg)	Max Speed (r/min)	MaxRCF ( xg)
12x1.5/2.2ml	16000	17800	16000	17800
18x1.5/2.2ml	----	----	13200	16110
24x1.5/2.2ml	13200	16110	13200	16110
36x0.5ml	16000	23669	13200	16110
10x5ml	13000	11400	13000	11400
12x10ml	12000	14800		
6x50ml	12000	13000		