



Product Name: Magnetic Stirrer with Hotplate

Model: MS300HS

Brand: MTOPS

Country of Origin: Korea

Features:
Various sizes of the top plate and flexible choice
Ceramic coated top plate provides excellent chemical-resistance, thermal-proof, and corrosion-resistance
Reflective, bright white color ceramic plate enhances sample visibility
Heating temperature is controlled by analog temperature controller(Built-in K type sensor)
High-heat conductivity of top plate provides thermal uniformity over the entire heating surface
Electronic feedback control maintains a precise speed
Smooth acceleration of stirring speed ensures no spillage

Specification:

Model		HS12	HS15	MS300HS	HS33
Type		Hot plate & Magnetic Stirrer			
Temperature	Range	Max 380°C			Max 350°C
	Display	Scale(temperature)			
	Controller	Electronic energy regulator			
Heater power		350W	480W	680W	1400W
Stirring	Speed	100~1500rpm			100~1300rpm
	Speed display	Scale(rpm)			
	Capacity(H ₂ O)	up to 1L	up to 2L	up to 5L	up to 15L
Speed Controllerer		Feedback control			
Top Plate		Ceramic coated top plate			
Dimensions	Plate area(mm)	120×120	150×150	180×180	300×300
	Overall(W×D×H)(mm)	160×280×107	160×280×107	200×310×107	310×430×120
Electrical supply		220V 50/60Hz or 110V 50/60Hz			

manager@psc.com.bd

www.psc.com.bd

Office address:

31, Manun Plaza (Ground, 1st & 2nd Floor)
Shahid Nazrul Islam Sharak, Tikatuli, Dhaka-1203

Weight	2.5Kg	2.6Kg	3.6Kg	6.8Kg
--------	-------	-------	-------	-------

Warranty: 01 (One) year warranty from the date of installation.

Installation & training: Free installation, operation & maintenance training will be provided at your laboratory site.

Paradise Sales Terms & Condition:

- **Delivery:** From ready stock / within 07 (seven) days. Without ready stock 6-8 weeks, after work order.
- **Payment:** By cash/Cheque

Contacts:

Sales : +88 01733-229381-01733-229386, 01733-229388
01713-034889

Service : +88 01733-229382, 01733-229391, 01733-229392
01733-229389

Accounts: +88 01733-229380

manager@psc.com.bd

www.psc.com.bd

Office address:

31, Manun Plaza (Ground, 1st & 2nd Floor)
Shahid Nazrul Islam Sharak, Tikatuli, Dhaka-1203