



SF6

Features

- ▶ For flocculation test of water, possible to stir 4 or 6 samples at the same time.
- ▶ The fixed illuminator provides glare-free lighting to provide diffused lighting of samples
- ▶ Life time keeps long due to the excellent & durable brushless DC motor.
- ▶ Speed control maintains 4 or 6 peddle's rotation speed concurrently by feedback control
- ▶ Digital speed indication & built-in timer(99hr 59min)
- ▶ Indicate separately set-in rpm an practical rpm on the bright LED screen and confirm at the same time.
- ▶ SP4 : portable type



SF4

Specifications

| Model | | SF6 | SF4 |
|--------------------|----------|------------------------|---------------|
| No. of places | | 6 | 4 |
| Stirring | Capacity | up to 2000ml | |
| | Motor | Brushless DC motor 60W | |
| | Speed | 20~300rpm | |
| Speed Controller | | Feedback control | |
| Speed Display | | Digital | |
| Timer | | 99hr 59min | |
| Dimensions (W×D×H) | | 900×250×450mm | 320×320×390mm |
| Electrical Supply | | 220V 50/60Hz | |
| Weight | | 21.5Kg | 9Kg |
| Price | | | |