



Features

- ▶ Homogenizer of rotor/stator type and possible to do quick Homogeneous and dispersion
 - ▶ Sample inserted into the rotor inside is circulated repeatedly through the narrow slit between rotor and stator, and emulsified, homogenized, mixed and pulverized thru dispersion by strong rotor.
 - ▶ Easy to change various rotor and stator by a driver only what you want to use
 - ▶ High quality tools dispersed durability (Stainless steel 316L)
 - ▶ easy to stick and take off for cleaning
- *Fixing clamp (option)

Specifications

Model	SR30
Motor(input/output)	670/340W
Working Capacity	10~2000ml
Speed Range	3700~27000rpm
Speed Adjustment	Stepless
Speed Display	Scale
Dimensions (W×D×H)	75×75×255mm
Weight	1.9Kg
Electrical Supply	AC 220V 50/60Hz
Price	
Price(set)	

※ Note : set ⇒ Homogenizer(SR30) + Tool(Ø20 Standard) + Stand(PS80) + Clamp(CH320)

Dispersing tools

Model	Shaft Ø	Price	
T8SF	Ø8mm		
T10SF	Ø10mm		
T20SC, T20CC, T20MC	Ø20mm		
T20SF, T20CF, T20MF, T20EF			
T30SC, T30CC, T30MC		Ø30mm	
T30SF, T30CF, T30MF, T30EF			
Rotor Ø20mm (Only)			
Rotor Ø30mm (Only)			