



- Powerful stirring and solid structure
- Large volume stirrer up to 50 liters
- Durable and powerful BLDC motor
- Max. speed : 1200rpm
- Digital display of speed with built-in timer
- Precise speed control with electronic feedback control
- Dual displays of the setting and actual rpm simultaneously
- Display of torque

### Features

- ▶ Powerful BLDC motor and "rare earth" magnet stirs large capacity stirrer
  - ▶ Smooth starting action provides to reach the desired speed smoothly
  - ▶ Digital display of speed with built-in timer(99hr 59min)
  - ▶ Electronic feedback control maintains a precise speed and torque
  - ▶ Clear LED display shows set and actual rpm
  - ▶ Easy to operate with touch key-pad
  - ▶ Stainless steel top plate is corrosive-resistance
  - ▶ Electronic safety circuit
- \*Supplied with dia. 12.7X76mm PTFE coated stirrer bar.

### Specifications

Model		MSM204D	MSM208D
Type		Magnetic stirrer	
Stirring	Capacity(H <sub>2</sub> O)	Max.50L	
	Motor	Brushless DC motor 20W	
	Speed	150~1500rpm	
	Speed display	Digital	
	Times	99hr 59min	
Speed Controller		Feedback control	
Plate		Stainless steel	
Dimensions	Plate area	360×360mm	
	Overall(W×D×H)	360×370×100mm	
Electrical Supply		220V 50/60Hz	
Weight		4.9Kg	
Price			