

PRECISE HEATING & STIRRING HEI-PLATE MIX 'N' HEAT MAGNETIC STIRRER

HIGH-RESISTANT HEATING PLATE

REFRACTORY HOUSING

800 W

MADE IN
GERMANY

3-year warranty on
all devices and an op-
erational lifespan of
more than 10 years

Hei-PLATE Mix 'n' Heat Core

Interface for more precision

For the EKT Hei-Con to measure and control the temperature precisely

Robust and durable for all challenges in the lab

Kera-Disk® plate – scratch-resistant, refractory, chemical-resistant and easy to clean

Innovative body design for additional protection against burns, secure handling and reliable trapping of liquids from the device

Refractory die-cast aluminum housing

Intuitive handling for a quick start

Standby button for switching the device easily on and off

Direct access to the main functions by two separate rotary controls

Heating button for the deliberate activation/deactivation of the heating function as well as for the visual display of the residual heat

Separate LED display of the heating



ACCESSORIES RECOMMENDATION



Silicone protective cover



Electronic contact thermometer



Heating and bath attachments



High-performance condenser