

KERN OCM 161

The inverted biological laboratory microscope

KERN



Category

Brand	Optics
Product category	Microscope
Product group	Inverted microscope
Product family	OCM-1

Construction

Dimension (W×D×H)	660×600×335 mm
Optical system	Infinity
Tube type	Trinocular
Tube angle of inclination	45°
Tube type	Siedentopf
Light distribution	0:100
Contrasting methods	Bright field Phase contrast Fluorescent (optional)

Standard objectives	10× 20× 40× PH 20×
---------------------	-----------------------------

Lens quality	Infinity Plan
Nosepiece screw-in locations	5
Interpupillary distance [Min]	48 mm
Interpupillary distance [Max]	76 mm
Diopter adjustment	both-sided
Diopter adjustment [Min]	-5
Diopter adjustment [Max]	5

Ocular

Ocular type	Eyepiece HWF 10x / Ø 22mm with anti-fungus, high eye point
Ocular field width	HWF

Eye point	High Eye Point
Ocular magnifications	10x
Ocular visual field	22 mm
Ocular diameter	30 mm

Objective

Objectives - Details	Objective Infinity Plan 10 x / 0,25 LWD, Fluor, anti-fungus Objective Infinity Plan 20 x / 0,45 LWD, Fluor, anti-fungus PH-Objective Infinity Plan 20 x / 0,45 LWD, Fluor, anti-fungus Objective Infinity Plan 40 x / 0,65 LWD, Fluor, anti-fungus
Objectives magnifications	10x / 20x / 40x

Illumination

Illuminance	Transmitted light
Illumination type transmitted light	Halogen
Illumination intensity transmitted light	30 W
Illumination dimmable	Transmitted light
Filter possible	✓
Brand of light bulb	Philips

Focussing

Field of view [Min]	0,55 mm
Field of view [Max]	2,2 mm
Torque regulation	✓
Fine drive minimum	0,002 mm
Focusing mechanism	coaxial coarse and fine drive

Condensor

Condensor type	Standard
----------------	----------

Power Supply

Supplied power supply	Power supply unit
Plug-in power supply type	Table power supply (separable)
Plug-in power supply / adapter for countries - included with the delivery	EURO
Plug-in power supply / adapter for countries - optional	EURO UK

KERN OCM 161



The inverted biological laboratory microscope

Input voltage power supply / power [Max]	100 - 240 V 50 / 60 HZ 1,5
--	-------------------------------

Input voltage device / power [Max]	12 V, 4 A
------------------------------------	-----------

Power socket for mains adapter	Hollow plug, inside plus, Ø outside 5,5 mm, Ø inside 2,5 mm, long 9,5 mm
--------------------------------	---

Environmental conditions

Storage temperature [Min]	-5 °C
---------------------------	-------

Storage temperature [Max]	40 °C
---------------------------	-------

Approval

CE mark	✓
---------	---

Packing & Shipping

Delivery time	1 d
---------------	-----

Dimensions packaging (W×D×H)	660×590×325 mm
------------------------------	----------------

Shipping method	Parcel service
-----------------	----------------

Net weight approx.	15 kg
--------------------	-------

Gross weight approx.	17 kg
----------------------	-------

Shipping weight	25,3 kg
-----------------	---------

Product information

GTIN/EAN number	4045761489738
-----------------	---------------

REACH classification	no prohibited substances
----------------------	--------------------------

Pictograms

STANDARD

